

# QuikClot®

## 10x10™ & 5x5™

### Hemostatic Dressing

## THE BLEEDING STOPS HERE.

*QuikClot® 10x10™ Hemostatic Dressing is a sterile, white, non-woven hemostatic gauze measuring 10 cm by 10 cm that is impregnated with an inert mineral, kaolin, that is both safe and effective in promoting clotting and does not contain any animal or human proteins. QuikClot® 10x10™ Hemostatic Dressing is packaged in a sterile, peel-able, foil pouch and should be used in conjunction with direct manual pressure. Additionally, QuikClot® 10x10™ contains an x-ray identifiable strip for easy detection.*

*QuikClot® 5x5™ Hemostatic Dressing is a sterile, white, non-woven hemostatic gauze measuring 10 cm by 10 cm that is impregnated with an inert mineral, kaolin, that is both safe and effective in promoting clotting and does not contain any animal or human proteins. QuikClot® 5x5™ Hemostatic Dressing is packaged in a sterile, peel-able, foil pouch and should be used in conjunction with direct manual pressure.*

### EFFECTIVE:

- Works Fast
- Published Clinical studies have proven the effectiveness of QuikClot® hemostatic agents\*
- Chosen by the U.S. Military as their exclusive choice for extreme bleeding in wartime injuries
- Active hemostatic ingredient (kaolin) has been recognized for more than 50 years as an effective clotting agent

### SAFE:

- No known contraindications and not invasive
- No human or animal derived proteins - no thrombin, fibrinogen or shellfish by-products
- Kaolin has been used safely in cosmetics for many years

### EASY TO USE:

- Little or no training required due to intuitive gauze format
- Flexible and absorbent - conforms readily to the site
- Will not break down or fall apart under pressure

### COST EFFECTIVE:

- Safer, yet less expensive than thrombin or shellfish based products
- Stops bleeding rapidly, saving from investing in other more expensive methods and therapies


\*Journal of Trauma - 9/2009 & 2/2010  
Journal of Interventional Cardiology

Z-MEDICA

# QuikClot® 10x10™ & 5x5™ Hemostatic Dressing

**THE BLEEDING STOPS HERE.**

## Thrombelastography (TEG) Clot Time



Source: University of California, Santa Barbara (modified)

A variety of extensive, independent, pre-clinical testing has been done separately at both the US Army's Institute for Surgical Research (USAISR) and the Naval Medical Research Center (NMRC).

As a result, QuikClot® Combat Gauze™ hemostatic agent was chosen #1 by the Department of Defense to be carried by all branches of the US Armed Forces, as a first-line treatment for life-threatening hemorrhage.


## DIRECTIONS FOR USE - QuikClot® 10x10™ & QuikClot® 5x5™




1. Open package and remove QuikClot®.



2. Place dressing over wound and use it to apply pressure directly over bleeding source.



3. Continue to apply pressure for 3 minutes or until bleeding stops.



4. Wrap and tie bandage to maintain pressure.

• **Item # 271**

**QuikClot® 5x5™  
Hemostatic Dressing**

Available in — 10 units per box

• **Item # 272**

**QuikClot® 10x10™  
Hemostatic Dressing**

Available in — 10 units per box

**To place an order:**

**Call: +1-203-294-0000**

**or**

**Fax: +1-203-294-0688**



**Z-MEDICA**